

Product Specification — CO-RE Tip

P/N and Product description

- P/N 194053
- 4ml 1T Rack CO-RE Filter Tip, non-sterile, black, for Microlab STAR Line



Material

- Polypropylen
- Filter tip insert: Polyethylene
- conductive

Packaging & processing

- 4 tips/tray, 24 trays/case, 96 tips/case
- Color-code label on tip tray: green
- Barcode label contains: unique identification number, Lot No., production date, cavity no.
- Barcode type: Code 128 subset B or C
- Tips usable on CO-RE single channel and 5000µl pipetting channels



Certification / Purity grade

- ISO Class 8 Clean Room production, compliant with ISO 9001, GMP and FDA
- CE IVD conform



- DNase-/RNase-free
- Pyrogen-free (endotoxin-free)
- PCR inhibitor-free

Labeling

Each box of disposable tips provides a label which indicates a barcode, part number (REF) and Lot No.

Measurement

- Tip length: 118mm
- Inner diameter tip orifice: 1.05mm

HAMILTON 

Web: www.hamiltonrobotics.com

Email: infoservice@hamiltonrobotics.com

United Kingdom & Ireland

Tel: +44 (0)121-717-0199

China

Tel: +86-21-6164-6567

France

Tel: +33 (01) 69751616

Italy

Tel: +39-39-689-33-93

Denmark, Norway,
Sweden, Finland

Tel: +46 (0) 8 410 273 73

Germany, Switzerland,
Austria, Benelux

Tel: +49 (089) 552649-0