

## Serological Pipette 1 – 70 ml










- Smooth interiors prevent cell hang-up or clumping inside
- Excellent optical transparency
- Optimized form of the mouthpiece that preserves the rubber of your pipette aid
- Volume identification by color coding on pipette, single-wrapping, dispenser, box label and cotton plugs to some extent
- Convenient, stackable TPP dispenser box
- lint-free paper/plastic single wrapping

### Additional features to the Serological Pipette with Reservoir:

- The reservoir of the 25 ml and 50 ml pipette augments the volume considerably without adding proportionate length to the pipette.
- Stable, break-resistant highly low drip pipette tip (picture 1)
- Slender, short form enables ergonomic working conditions in the sterile hood
- Small outside diameter of the tube at the tip resulting in an excellent accessibility of these large pipettes even in small tissue culture flasks
- Optimized form of mouthpiece
- Fits in pipettor Turbo-Fix or all commercially available pipette aids with rubber inserts or the Stabifix from TPP

### Features of the standard Serological Pipette:

- Accurate volume 1 – 70 ml
- Highly visible, bidirectional graduation
- Larger volume by over-volume graduation, still shorter ergonomic design for easier maneuvering

Type	Product-No.	Volume ml	Reservoir ml	Colour Code	Dimensions Ø x l mm	Volume Graduations ml	Material	Qty / Case
	94001	1		yellow	4.8 x 270	1/100	PS	400
	94002	2		green	5.5 x 270	1/100	PS	300
	94005	5		blue	9.5 x 295	1/10	PS	200
	94010	10		orange	11.0 x 295	1/10	PS	200
	94024	25		red	16.0 x 300	2/10	PS	100
	94525	25	15	red	15 x 345	2/10	PS	60
	94550	50	20	purple	20 x 345	1/10	PS	50

1 Small pipette tip.

2 TPP product line of serological pipettes.

