

**Bellco Glass, Inc.**  
**Specification Sheet**  
**Product #: 1985-85000**

**Description: Collector™ Complete, 85mL**

- **Tissue Sieve Homogenizer, used to prepare cell suspensions and gel extractions of proteins**
- **All 304 stainless steel construction**
- **Pan Dimensions: 1-1/4" OD x 2-1/2" W**
- **Complete consists of: Collector™ pan & key, glass pestle, 9 screens; 1 Ea from 10 to 150 mesh**



**The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.**



**BELLCO GLASS, INC.**



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05