



RNase-free water

Cat. No.	size
E0210-01	5 x 1.5 ml
E0210-02	5 x 5 ml
E0210-03	500 ml

Storage Conditions: Store at 2-8°C

DEPC-treated and nuclease-free water for RNA work.

Description:

EURx Diethylpyrocarbonate (DEPC)-treated water is guaranteed nuclease-free (tested for both DNase and RNase) and is suitable for all RNA work.

Quality Control:

Each lot is assayed for nuclease activity.

This product is developed, designed and sold exclusively for research purposes and in vitro use only.