

# Technical Data Sheet

Doc No – Ab/TDS/GLP-157-I

Rev No : 07 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**PP- Beaker without Handle**

**Valid for Catalog No. P50601, P50602, P50603, P50604, P50605, P50606, P50607, P50608**

### **1. Description/Specification**

#### **1.1 Description:**

Reusable Polypropylene (PP) Beaker without Handle

1.1.1 Tapered stacking style saves space and makes it easy to handle.

1.1.2 All sizes have easy to read silk-screened graduations.

1.1.3 Precise pouring spouts.

#### **1.2 Design:**

Tapered Stacking Design

#### **1.3 Performance Specifications:**

##### 1.3.1

Code	Capacity (ml)	Printing (Min - Mix)	Sub division	Height (mm)	No./ Case
P50601	50 ml	10ml - 50ml	5ml	59 mm	12
P50602	100 ml	20ml - 100ml	10ml	73 mm	12
P50603	250 ml	50ml -250ml	25ml	94 mm	12
P50604	500 ml	100ml - 600ml	50ml	128 mm	6
P50605	1000 ml	100ml - 1000ml	50ml	150 mm	4
P50606	2000 ml	200ml - 2000ml	100ml	192 mm	4
P50607	5000 ml	1000ml - 5000ml	500ml	263 mm	2
P50608	10000 ml	2000ml - 10000ml	1000ml	328 mm	1

1.3.2 Tolerance limit  $\pm 10\%$

#### **1.4 Material:**

1.4.1 PP- Polypropylene (PP), confirm US FDA 21CFR

1.4.2 Complies with the requirement of European Pharmacopoeia, Plastic Containers and closures for Pharmaceutical use.

1.4.3 The product has passed the USP testing including class VI tests ,Materials has successfully passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email: labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**

# Technical Data Sheet

1.4.4 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations 1935/2004/CE, 2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & its following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water

1.4.5 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

## 1.5 Color:

Translucent (Natural)

## 1.6 Graduations:

1.6.1 Large, Printed with temperature, chemical and solvent resistant ink.

1.6.2 Meets accuracy requirements under "DIN 12681/ISO 7056 – 1981 (E) & BS 5404 Part 1".

## 1.7 Sterilization:

NA

## 1.8 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

## 1.9 Other Information:

Reusable, can be autoclaved up to 10 autoclave cycle.

## 2. Basic Features

2.1 Made up of High purity Virgin Polypropylene (PP)

2.1 Temperature range:  $-80^{\circ}\text{C}$  to  $+121^{\circ}\text{C}$

2.3 Autoclavability: Yes, Autoclavable at 15 psi and  $121^{\circ}\text{C}$  for 20 minutes

2.4 Chemical resistance: See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

2.5 Manufactured in Aseptic Environment

## Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at  $121^{\circ}\text{C}$  for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least  $50^{\circ}\text{C}$  before removing product.

# Technical Data Sheet

## Instruction for Use

Product can be used for food contact also.  
Storage temp. Max 30°C  
Before using the product check crackness etc.  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimen.

## WARNING

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure  
Consult local fire codes prior to storage of flammable liquids.

**UKRP**  
Obelis UK LTD  
Sandford Gate,  
East Point Business Park,  
Oxford OX4 6LB  
United Kingdom  
Tel: +44 1891 378012  
info@obelis.co.uk

**Obelis s.a.**  
Bd General Wahis 53  
1030 Brussels Belgium  
Tel: +32 2 732-59-54  
Fax: +32 2 732-60-93  
mail@obelis.net

**NON STERILE**

30°C

**IC REP.**

**Abdos Labtech Pvt. Ltd.**  
Plot No. 1 Lakshewari Bhagwanpur Roorkee 247667 (UK) India  
Email : labtech@abdosindia.com, www.abdoslifesciences.com  
Ph: +(91) 11 41081897

Bictronic IFU

