



CHEMWARE® GAS BAGS & LIQUID SAMPLING

FEP FLUOROPOLYMER FILM

- Chemically resistant
- Flexible
- Will not contaminate sample

Made from 5 mil thick FEP fluoropolymer film. These flexible, sturdy bags are heat-sealed on all four sides for complete chemical resistance and leakproof operation. Inert to almost all chemicals. Low permeability of FEP limits sample loss. Available with either a nickel-plated, brass on/off valve that accepts 1/4-inch I.D. flexible tubing or with PTFE/silicone septum in polypropylene housing for use with a hypodermic needle.

Easy to use. Unlimited storage life prior to use. Can be used for gas sample analysis, extracting air samples from smokestack exhaust systems, monitoring air quality or calibrating instruments in the field. Individually packaged. Custom sizes and configurations available on request.

Catalog No.	Fitting	Size in.	Capacity L
D1075002	On/Off	6x6	0.5
D1075001	On/Off	6x10	1.2
D1075003	On/Off	12x12	4.7
D1075004	On/Off	24x36	50.0
D1075005	Septum	6x6	0.5
D1075006	Septum	6x10	1.2
D1075007	Septum	12x12	4.7
D1075008	Septum	24x36	50.0

All bags also available in case of 10, except 50L size which is available in a case of 5.

CHEMWARE® GAS BAGS & LIQUID SAMPLING

PVF FILM

- EPA cited
- Flexible
- Will not contaminate sample
- Choice of fitting style: on/off or septum
- For use in gas sampling trains
- Use for gas or liquid sample containment

Reusable bags made of 2 mil thick PVF fluoropolymer film. An economical alternative for analysis of non-corrosive gases. Durable, heat-sealed seams ensure leakproof operation. Available with either a nickel-plated, brass on/off valve that accepts 1/4-inch I.D. flexible tubing or with a PTFE/silicone septum in a polypropylene housing for use with a hypodermic needle.

Uses include sampling carbon monoxide, sulfur dioxide, hydrogen sulfide, radon and mercaptans. Cited in Title 40 CFR, PT.60, App A Meth. 18. Individually packaged. When used for liquid samples, bags should be supported with a tray. Not pressure rated. Custom sizes and configurations available on request.



LATEX GAS BAGS & LIQUID SAMPLING

LATEX NATURAL RUBBER

- Pure rubber
- Various capacities

Cylindrical with tube and ebonite valve.

Catalog No.*	Capacity	Pack Size
101912	0.5 liter	10
101914	1 liter	10
101916	2 liters	10
101920	5 liters	5
101922	10 liters	5

20, 30 and 50 liters also available.

Bags size 0.3 to 8.1, also available in case of 10.
Bags size 20.3 to 85.7, also available in case of 5.

Catalog No.	Fitting	Size in.	Capacity L
D1075009	On/Off	6x6	0.3
D1075010	On/Off	9x9	1.6
D1075011	On/Off	12x12	3.8
D1075012	On/Off	15x15	8.1
D1075013	On/Off	18x18	20.3
D1075014	On/Off	24x24	37.7
D1075015	On/Off	30x30	85.7
D1075016	Septum	6x6	0.3
D1075018	Septum	9x9	1.6
D1075017	Septum	12x12	3.8
D1075019	Septum	15x15	8.1
D1075020	Septum	18x18	20.3
D1075021	Septum	24x24	37.7
D1075022	Septum	30x30	85.7