

Pyrex® Evaporating Dish, with spout, 90 mL

- Complies to DIN 12336
- Manufactured from borosilicate glass 3.3
- Curved sides with spout
- Flat base
- High resistance to chemical attack and thermal shock
- 1450/01D, 1450/03D & 1450/04D are unbadged.

Catalog Number: **1450/04D**

Product Overview

Transferring substances and liquids into other states of aggregation and thus initiating further reactions is very important for the daily work in the laboratory. A wide variety of vessels, apparatus and equipment are used for these work.

Specifications

Capacity (mL)	90
Diameter (mm)	80
Height (mm)	45
Quantity	10

This information is accurate as of 16/07/2024 and is subject to change.