

PYREX[®] Measuring Cylinder, Class A, UKAS certified, 500 mL

- Complies with ISO 4788
- Manufactured from borosilicate glass 3.3
- Calibrated to Class A tolerances
- Individually tested and calibrated in accordance with the requirements of ISO 17025 and accredited by the United Kingdom Accreditation Service (UKAS)
- Supplied with a UKAS calibration certificate showing the exact volume contained at 5 points on the graduated scale
- Conformity certified
- Inscribed with an individual serial number for traceability purposes
- Graduations in white enamel

Catalog Number: **3201/14UKAS**

Product Overview

These amber PYREX[®] volumetric flasks are tested and calibrated to the highest standards in accordance with the requirements of the United Kingdom Accreditation Service (UKAS).

Specifications

Tolerance (mL)	± 2.50
Capacity (mL)	500
Diameter (mm)	53.35
Height (mm)	365
Graduation Intervals (mL)	5
Quantity	1

