

Specifications

Unit Dimensions(WxLxH)	18.5 x 22 x 3cm
Gel Viewing Dimensions(WxL)	12 x 18cm
Input Voltage	100-240Vac
Input Current	2.0A
LED Source	Built-in blue light & white light LED module
LED life(hours)	>30,000
Emission maxima(nm)	470nm
Store temperature	25°C
Automatic Power-Off	5 mins
Filter Type	Amber filter (580nm)
Certifications	CE/cUL

Related Ordering Information

Cat. No.	Description
BP001CU	BluPAD Dual LED Blue/White Light Transilluminator
BP-miniDR	BluPAD miniDarkroom
BP-KNIFE	BluPAD Gel-Cutting Knife
BP-AFC	BluPAD Upper amber filter cover
BP-WUP	BluPAD white uniform plate
BP-BDPG	BluPAD Blue uniform plate w glass
BP-AC	BluPAD 60 watt adaptor

Manufacturer

BIO-HELIX
CO., LTD.

Bio-Helix Co., Ltd.

Room MAF 312, No.2, Beining Rd., Keelung City Taiwan (202), TAIWAN

Tel: +886-2-24624956


E-mail: info@bio-helix.com



www.bio-helix.com

10th 2007-2017
ANNIVERSARY

B/uPAD



Dual LED Blue/White Light Transilluminator