


**DATASHEET****PART NUMBER****257305****SolidFrame PCR plates****Description :**

- The ultra Rigid fully skirted frame remains dimensionally stable during thermal cycling, which facilitates reliable robotic handling and ensures that the seal to each well is not compromised during thermal cycling
- Thin-walled tubes made from 100% virgin medical grade polypropylene
- Rigid polycarbonate frame
- Low profile wells have a maximum capacity of 225 µl
- Manufactured to the highest quality standards for your PCR and qPCR applications
- Our fully skirted 2-component plate fits a wide range of instruments and is fully automation compatible
- Certified free from Dnase, Rnase, nucleases and human gDNA

**Specification :**

- Number of well : 96
- Frame color : black
- Well color : white high sensitive
- Skirt : skirted
- Style : standard
- Well volume : 0,1 ml

**Units : 50 plates**