

ClearLine membranes filters

The MCE membranes are made with cellulose acetate and cellulose nitrate fibres, they are naturally hydrophilic, mechanically stable and have a high loading capacity. They are suitable for microbiological analysis, for colonies counting or for pre-filtration and clarification of samples.

Membrane : MCE - Mixed cellulose ester

Pore size (for Ø 47 mm) : 0,45 µm

Water flow @ 10 psi : 38-71 ml/min/cm²

Bubble point : 29-40 psi

Diameter : 47 mm

Sterile

Colour : Black, white gridded

Packaging : 4 bags of 25

**Units :**

100 membranes