

## Certificate of EU Medical Device Notification

This is to certify that, according REGULATION (EU) 2017/746 of THE EUROPEAN PARLIAMENT AND OF THE COUNCIL on in vitro diagnostic medical devices,

**HUMISS TRADING S. L.**

**Address: Calle Luis Buñuel 12, MADRID, 28018, Spain**

**SRN: ES-AR-000040177**

has fulfilled all notification responsibility and duty as the European Authorized representative of:

**Manufacturer: Suzhou Yuanzan Intelligent Manufacturing Technology Co., Ltd.**

**Address: Building A1/A5, No. 188 East Shengxiang Road, Qiandeng Town, Kunshan City, Jiangsu Province, 215341, P.R. China**


The manufacturer has provided with all the appropriate declaration according the REGULATION (EU) 2017/746 requirements including the EC Declaration of Conformity confirming that the in vitro diagnostic medical devices, as stipulated here below, is fulfilling the applicable requirements of REGULATION (EU) 2017/746.

**Product(s) details shown in the Annex of this certificate.**

Where then manufacturer affixes the CE marking to the product listed, they must ensure that all the requirements of the appropriate EU regulation(s) have and continue to be met.

The notification of aforementioned device(s) has been completed by European Authorized representative in Spain, the Spain Competent Authority has notified the manufacturer's in vitro diagnostic medical device above and has allocated registration.

You can check the authenticity of the document at the headquarters of the AEMPS (<https://localizador.aemps.es>) with CSV number listed above.



*This is only a CE mark sample which is only use for reference.*



Signature of Executive Director

Tamara Manchón

<b>Product name</b>	<b>Models</b>	<b>Activities</b>	<b>Registration number</b>	<b>CSV</b>	<b>Classification</b>
Pipette Tips	5µl, 10µl, 20µl, 30µl, 50µl, 70µl, 100µl, 200µl, 250µl, 300µl, 500µl, 700µl, 1000µl, 1250µl, 2000µl, 5000µl	Manufacturer	RPS/99/2026	BZALGBZB96 SR9DHP956D	Class A, Rule 5
Centrifuge Tube	15ml, 50ml, 250ml, 500ml	Manufacturer			Class A, Rule 5
Microcentrifuge Tubes	0.6ml, 1.5ml, 2.0ml, 5.0ml	Manufacturer			Class A, Rule 5
Cryogenic Vials	0.5ml, 1ml, 1.2ml, 1.5ml, 2.0ml, 3.0ml, 4.0ml, 5.0ml	Manufacturer			Class A, Rule 5
Cell Culture Flasks	12.5cm <sup>2</sup> , 25cm <sup>2</sup> , 75cm <sup>2</sup> , 175cm <sup>2</sup> , 182cm <sup>2</sup> , 225cm <sup>2</sup> , 300cm <sup>2</sup> , 870cm <sup>2</sup>	Manufacturer			Class A, Rule 5
Culture Dishes	35mm, 60mm, 65mm, 70mm, 90mm, 100mm, 150mm	Manufacturer			Class A, Rule 5
Cell Culture Plates	1, 4, 6, 12, 24, 48, 96, 8-strip, 384	Manufacturer			Class A, Rule 5
Serological Pipettes	2ml, 5ml, 10ml, 25ml, 50ml, 100ml, 0.5ml, 3ml, 4ml, 6ml, 7ml, 8ml, 9ml, 15ml, 20ml, 30ml, 40ml	Manufacturer			Class A, Rule 5
Deep-Well Plates	0.5ml, 1.0ml, 1.6ml, 2.2ml, 10ml	Manufacturer			Class A, Rule 5
Tip Combs	24, 48, 96, 8-Strip	Manufacturer			Class A, Rule 5
PCR Tubes	0.1ml, 0.2 ml	Manufacturer			Class A, Rule 5
PCR 8-Strip Tubes	0.1ml, 0.2 ml	Manufacturer			Class A, Rule 5
PCR Plates	40ul, 0.1 ml, 0.2 ml	Manufacturer			Class A, Rule 5
Elisa Plates	8, 12, 16, 48, 96	Manufacturer			Class A, Rule 5
Pipettes	Single channel pipettes; 8-channel pipettes; 12-channel pipettes	Manufacturer			Class A, Rule 5
<i>End of the annex.</i>					