


1000 µl tip, LTS compatible, graduated, sterile, MultiRack, blue tray
ClearLine

- Volume : 1000 µl
- Longueur : 98,6 mm
- Packaging : In rack Clearline Plus - 8 x 96
- Compatibility : Compatible with Rainin LTS™ micropipettes
- Tray color : blue
- Low retention : no
- Autoclavable : yes
- Sterile : yes
- Filter : no
- Color : natural
- Graduation : 200, 500, 1000
- Material : polypropylene



Units : 768