


### Shorty serological pipettes

- Manufactured from highly transparent polystyrene (PS) material
- With a filter plug for filtration
- Different color bands
- Clear dual-scale lines for positive and negative readings
- Easy to handle under a cabinet
- Same precision as a standard pipette
- Free from DNA enzymes, RNA enzymes, and pyrogene
- Irradiation sterilization
- Package : Bulk
- Volume : 10 ml
- Graduation : 0,200000000000000001
- Color code : Orange
- Length : 232 mm
- Packing : 2 bags of 100

**Units :**

200