

**Bellco Glass, Inc.  
Specification Sheet**

**Product #: 5630-12000**

**Description: Roux Culture Bottle with Delong Style Neck**

**Material: Duran® Borosilicate Glass type 3.3 as specified in international standard ISO 3585.**

**Capacity: 1200mL**


**Overall Dimensions (WxDxH): 121x57x290mm**

**Neck Size: 38mm OD Delong Style designed for use with 38mm closures**

**Offset Neck Style**

**Application/Usage: This bottle is designed with flat surface bottom when laid down horizontally holds approximately one-half total volume. Dam in neck prevents culture from entering neck area. It is useful for growing mass cultures and single or monolayer cultures.**

**Sterilization: Yes, wet or dry**



**The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.**



**BELLCO GLASS, INC.**



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05