

**Bellco Glass, Inc.
Specification Sheet**

Product #: 5632-12000

Description: Roux Culture Bottle with Screw Neck

Material: Duran® Borosilicate Glass type 3.3 as specified in international standard ISO 3585.

Capacity: 1200mL

Overall Dimensions (WxDxH): 132x56x280mm

Neck Size: 38mm Screw Thread

Screw Cap Included

Application/Usage: This bottle is designed with flat surface bottom when laid down horizontally holds approximately one-half total volume. Dam in neck prevents culture from entering neck area. It is useful for growing mass cultures and single or monolayer cultures.

Sterilization: Yes, wet or dry



The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.



BELLCO GLASS, INC.



Bellco Glass, Inc. · 340 Edrudo Road · Vineland, NJ · 08360

Ver: 7/05