

# Technical Datasheet

---

**Article No. 2036**

**Glycerol anhydrous p. a.**  
**min. 99.5 % C<sub>3</sub>H<sub>8</sub>O<sub>3</sub>**

**For laboratory use.**

Parameter	Value
Appearance	Colourless liquid
Density	1.26 g/ml
Mol.-Weight	92.09 g/mol
Melting point	18 °C
Boiling point	290 °C
Water	max. 0.5 %
Colour (Hazen)	max. 8
Sulfated ash	max. 0.01 %
Heavy metals (as Pb)	max. 0.0005 %
Chloride (Cl)	max. 0.001 %
Hylogen compounds	max. 0.035 %