

## Product Specification

<b>Catalogue No:</b>	1130/02M
<b>Product Description:</b>	50ml Conical Flask, Narrow Neck
<b>Brand:</b>	PYREX
<b>Type of Glass:</b>	Borosilicate glass Type 3.3, which conforms to ISO3585 European and United States Pharmacopoeia Type 1 glass
<b>Product Markings:</b>	Capacity figures, graduation lines, product branding, marking spot and a trace-code (Colour white) A batch certificate is available from our website: <a href="http://www.scilabware.com">www.scilabware.com</a>
<b>Dimensions:</b>	Overall height = 85mm approx. External diameter = 51mm approx. Internal diameter of the neck = 18mm approx.
<b>Product Standard:</b>	The product complies with the requirements of ISO1773 and DIN 12353
<b>Pack Quantity:</b>	Quantity 10 pieces per pack

SciLabware Limited confirms that the above mentioned product has been manufactured to meet all relevant standards as detailed by our Quality Management System.



**R J Turner**  
(Quality Manager)