



Large volume container - 100 L - white

Ref. : G140521



Plus points

Large volume container.
Excellent mechanical resistance: made from low density polyethylene for increased shock resistance.
The durable, thick walls mean no distortion and warping.
Limited gaps and overlapping, equipped with edges so easy to handle.

Product features

Height	485 mm
Width	535 mm
Length	690 mm
Capacity	100 L
Material	PEHD
Colour	White
Product bottom internal length	550 mm
Product bottom internal width	400 mm
Product bottom internal height	470 mm
Weight	6.300 kg
GENCODE	3100411405214
Customs code	39231000

Logistic data

Number of items per pallet	15
Pallet weight	104 kg
Pallet length	800 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.