



Barrel + metal toggle joint closure - 100 L - white

Ref. : G178921



Plus points

Food contact certified plastic barrel.
Suitable for a large number of uses (brining, pickling etc.) when transporting and storing liquid or solid food products.
Global and specific migration test approved.
Contents are protected from contamination thanks to its watertight lid. .
Handling is easy and can be opened completely.

Product features

Height	680 mm
Diameter	492 mm
Capacity	100 L
Material	PEHD
Colour	White
Product internal diameter	400 mm
Weight	5.140 kg
GENCODE	3100411789215
Customs code	39233090

Logistic data

Number of items per pallet	12
Pallet weight	71 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	2070 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.