



Lid for 8 L flat tray - clear

Ref. : G160341



Plus points

Dual function lid for flat trays.
Closing & ventilation preventing odours from escaping.
Contents can be seen through the lid.

Product features

Height	20 mm
Width	340 mm
Length	490 mm
Material	PP
Colour	Transparent
Weight	0.330 kg
GENCODE	3573671603414
Customs code	39235090

Logistic data

Number of items per pallet	400
Pallet weight	186 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.