



Lid for 3 L flat tray - clear

Ref. : G160151



Plus points

Dual function lid for flat trays.
Closing & ventilation preventing odours from escaping.
Contents can be seen through the lid.

Product features

Height	15 mm
Width	240 mm
Length	355 mm
Material	PP
Colour	Transparent
Weight	0.150 kg
GENCODE	3573671601519
Customs code	39235090

Logistic data

Number of items per pallet	1584
Pallet weight	282 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.