



## Lid for 25 L 35 L and 55 L containers - grey

Ref. : G193377



### Plus points

Lid for rectangular plastic food container.  
Attached lid for increased hygiene.

### Product features

Height	35 mm
Width	400 mm
Length	600 mm
Material	PEHD
Colour	Grey
Weight	0.860 kg
GENCODE	3100411933779
Customs code	39235090

### Logistic data

Number of items per pallet	220
Pallet weight	189,2 kg
Pallet length	800 mm
Pallet width	1200 mm
Pallet height	1500 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.