




Rectangular container - 15 L - white

Ref. : G119621



Plus points

400x300 rectangular plastic food container

For sub-zero and ambient temperature storage and transporting raw materials.

No gaps nor overlaps.

Machine washable.

Attached lid for increased hygiene, optional extra.

Stackable without lid. Compatible with all GILAC 4 and 6 wheel dollies.

Product features

Height	235 mm
Width	300 mm
Length	400 mm
Capacity	15 L
Material	PEHD
Colour	White
Product bottom internal length	335 mm
Product bottom internal width	235 mm
Product bottom internal height	210 mm
Weight	1.090 kg
GENCODE	3100411196211
Customs code	39231000

Logistic data

Number of items per pallet	120
Pallet weight	142 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1850 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.