




Rectangular container - 55 L - white

Ref. : G119921



Plus points

600x400 rectangular plastic food container.
For sub-zero and ambient temperature storage and transporting raw materials.
No gaps nor overlaps.
Machine washable.
Attached lid for increased hygiene, optional extra.
Stackable without lid. Compatible with all GILAC 4 and 6 wheel dollies

Product features

Height	315 mm
Width	400 mm
Length	600 mm
Capacity	55 L
Material	PEHD
Colour	White
Product bottom internal length	520 mm
Product bottom internal width	330 mm
Product bottom internal height	315 mm
Weight	2.200 kg
GENCODE	3100411199212
Customs code	39231000

Logistic data

Number of pieces by subpacking unit	45
Subpacking weight	97,2 kg
Number of items per pallet	45
Pallet weight	108 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	2010 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.