



Lid for 12 L and 15 L containers - white

Ref. : G118921



Plus points

Lid for rectangular plastic food container.
Attached lid for increased hygiene.

Product features

Height	30 mm
Width	300 mm
Length	400 mm
Material	PEHD
Colour	White
Weight	0.320 kg
GENCODE	3100411189213
Customs code	39235090

Logistic data

Number of items per pallet	384
Pallet weight	198,2 kg
Pallet length	800 mm
Pallet width	1200 mm
Pallet height	1700 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.