



Grid for 8 L tray

Ref. : G100221



Plus points

Grids for flat trays.

Highly durable and shock resistant.

Enables the use of flat dishes as defrosting trays by separating the ice water from the product being defrosted.

Product features

| | |
|--------------|---------------|
| Height | 8 mm |
| Width | 270 mm |
| Length | 420 mm |
| Material | PEHD |
| Colour | White |
| Weight | 0.170 kg |
| GENCODE | 3100411002215 |
| Customs code | 39241000 |

Logistic data

| | |
|----------------------------|---------|
| Number of items per pallet | 1760 |
| Pallet weight | 346 kg |
| Pallet length | 1000 mm |
| Pallet width | 1200 mm |

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.