




## HACCP flat tray + lid + grid - 8 L - red

Ref. : G190334



### Plus points

Flat tray Complete package. Tray + Shelf+ Lid  
Dimensions suitable for refrigerators  
Can be used as :

- defrosting trays thanks to its wide dividing grid that separates the ice water from the product being defrosted .
  - preparation container for products that drip ingredients tray,
  - containers for preserving comparison samples,
  - Fermenting and decanting containers.
- Nesting and stacking with lid.  
Traceability management thanks (IML) injection moulding.

### Product features

Height	90 mm
Width	337 mm
Length	491 mm
Capacity	8 L
Material	PEHD
Colour	Red
Product bottom internal length	420 mm
Product bottom internal width	265 mm
Weight	1,130 kg
GENCODE	3573678674219
Customs code	39241000

### Logistic data

Number of items per pallet	98
Pallet weight	120,74 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1380 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.