



Grid for 10 L tray

Ref. : G100321



Plus points

Grids for flat trays.
Highly durable and shock resistant.
Enables the use of flat dishes as defrosting trays by separating the ice water from the product being defrosted.

Product features

Height	8 mm
Width	320 mm
Length	470 mm
Material	PEHD
Colour	White
Weight	0.210 kg
GENCODE	3100411003212
Customs code	39241000

Logistic data

Number of items per pallet	1260
Pallet weight	309 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.