



Food specimen containers + lids - 3 L - Pack of 5

Ref. : G160450



Plus points

"Pack of ""food specimen"" containers.
Containers especially for storing food specimen samples complying with HACCP catering standards.
Markings do not wash off, even in the machine."

Product features

Height	70 mm
Width	235 mm
Length	350 mm
Capacity	3 L
Material	PEHD
Product bottom internal length	280 mm
Product bottom internal width	165 mm
Weight	2.200 kg
GENCODE	3573678665903
Customs code	39241000

Logistic data

Number of items per pallet	143
Pallet weight	314,6 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.