


## GILACTIV antibacterial container + lid - 12 L - blue

Ref. : G119690



### Plus points

GILACTIV 400x300 antibacterial plastic food container.  
Proven antibacterial efficiency on the container surface against a wide range of germs, full or empty.  
No gaps nor overlaps and no sharp corners.  
Rewritable identification cartridge for inventory tracking.  
Recessed handles for better handling.  
Smooth interior walls and rounded corners for easier cleaning.  
Machine washable.  
Compatible with all GILAC 4 and 6 wheel dollies.

### Product features

Height	235 mm
Width	300 mm
Length	400 mm
Capacity	15 L
Material	PEHD
Colour	Blue
Product bottom internal length	335 mm
Product bottom internal width	235 mm
Product bottom internal height	210 mm
Weight	1.510 kg
GENCODE	3573678669000
Customs code	39231000

### Logistic data

Number of items per pallet	80
Pallet weight	130 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1890 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.