


## GILACTIV antibacterial container + lid - 35 L - blue

Ref. : G119890



### Plus points

- GILACTIV 600x400 antibacterial plastic food container.
- Proven antibacterial efficiency on the container surface against a wide range of germs when full or empty.
- No gaps nor overlaps and no sharp corners.
- Rewritable identification cartridge for inventory tracking.
- Recessed handles for better handling.
- Smooth interior walls and rounded corners for easier cleaning.
- Machine washable.
- Compatible with all GILAC 4 and 6 wheel dollies.

### Product features

Height	230 mm
Width	400 mm
Length	600 mm
Capacity	35 L
Material	PEHD
Colour	Blue
Product bottom internal length	530 mm
Product bottom internal width	335 mm
Product bottom internal height	210 mm
Weight	2.460 kg
GENCODE	3573678669024
Customs code	39231000

### Logistic data

Number of items per pallet	45
Pallet weight	121 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	2105 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.