




## Solid dough container - 15 L - white

Ref. : G179421



### Plus points

600x400 solid crate.  
For caterers and restaurateurs: storage crate for raw materials.  
Stacking and without lid can be stacked to a great height  
Easy cleaning flat bottom.

### Product features

Height	90 mm
Width	400 mm
Length	600 mm
Capacity	15 L
Material	PEHD
Colour	White
Product bottom internal length	565 mm
Product bottom internal width	365 mm
Product bottom internal height	80 mm
Weight	1,120 kg
GENCODE	3100411794219
Customs code	39231000

### Logistic data

Number of items per pallet	100
Pallet weight	140 kg
Pallet length	800 mm
Pallet width	1200 mm
Pallet height	2130 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.