




## Round container - 1519 - 50 L - white

Ref. : G030621



### Plus points

Round food container

### Product features

Height	525 mm
Diameter	450 mm
Capacity	50 L
Material	PEHD
Colour	White
Weight	3.090 kg
GENCODE	3100410306215
Customs code	39231000

### Logistic data

Number of items per pallet	60
Pallet weight	194,4 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.