



## Lid for 50 L container 1519 - white

Ref. : G152021



### Plus points

Lid for round G030621 food containers (G030621) and for 35 L small tubs (G608121).

### Product features

Height	35 mm
Diameter	495 mm
Material	PEHD
Colour	White
Weight	0.700 kg
GENCODE	3100411520214
Customs code	39235090

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.

### Logistic data

Number of pieces by subpacking unit	23
Subpacking weight	18 kg
Subpacking length	600 mm
Subpacking width	500 mm
Subpacking height	400 mm
Number of items per pallet	368
Pallet weight	298 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1700 mm
Number of units per layer	4
Number of layers per pallet	4