



Jerrycan - 20 L - white

Ref. : G105421



Plus points

Jerrycan fit for food contact.
170 mm spout.
Can be carried using the handle.
Highly durable.
Translucent meaning that the liquid level can be monitored.
Not suitable for storing hydrocarbons.

Product features

Height	460 mm
Width	200 mm
Length	320 mm
Capacity	20 L
Material	PEHD
Colour	White
Product bottom internal length	300 mm
Product bottom internal width	190 mm
Weight	1.435 kg
GENCODE	3573678666610
Customs code	39233090

Logistic data

Number of items per pallet	64
Pallet weight	100 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1970 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.