



Tote caddy - yellow

Ref. : G161428



Plus points

Tote caddy. Highly durable.
Suitable for universal use.
Tote caddy can be used with 1 or 2 compartments.

Product features

Height	205 mm
Width	255 mm
Length	385 mm
Material	PEHD
Colour	Yellow
Product bottom internal length	360 mm
Product bottom internal width	235 mm
Product bottom internal height	110 mm
Weight	0.840 kg
GENCODE	3573678670280
Customs code	39231000

Logistic data

Number of items per pallet	112
Pallet weight	116 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1500 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.