



## Tote caddy divider - Yellow

Ref. : G161228



### Plus points

Tote caddy divider.

### Product features

Height	42 mm
Width	104 mm
Length	382 mm
Material	PEHD
Colour	Yellow
Weight	0.200 kg
GENCODE	3573678670297
Customs code	39231000

### Logistic data

Number of items per pallet	1120
Pallet weight	164 kg
Pallet length	1000 mm
Pallet width	1200 mm
Pallet height	1460 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.