



Square dustbin - 4 L - white

Ref. : G038421



Plus points

Dustbin suitable for institutional catering and hotel bedroom use etc.
With interior white plastic bucket + handle
Highly durable
Stable and easy to maintain.

Product features

Height	228 mm
Width	230 mm
Length	200 mm
Capacity	4 L
Material	PP
Colour	White
Weight	0,470 kg
GENCODE	3573670384215
Customs code	39229000

Logistic data

Number of pieces by subpacking unit	144
Subpacking weight	67,68 kg

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.