



Carboy + metal toggle joint closure - 20 L - white

Ref. : G174721



Plus points

Carboy with screwtop. Migration tested and food contact certified guaranteeing quality product preservation.

Product features

Height	454 mm
Width	280 mm
Length	280 mm
Diameter	215 mm
Capacity	20 L
Material	PEHD
Colour	White
Weight	1.660 kg
GENCODE	3573678672116
Customs code	39233090

Logistic data

Number of items per pallet	60
Pallet weight	109 kg
Pallet length	1000 mm
Pallet width	1200 mm

All the technical information on temperature use is given as an indication for the products concerned. GILAC products are used at the user's own risk.