

Lowry Assay Kit

Pack Sizes

Certificates of Analysis (Lot.-no. - release date)

PDF Documents

Fast protein quantification with a sensitivity of 50 µg protein/ml, based on Lowry's method (1). It contains ready-to-use reagents and protein standard. Sufficient for 250 2-ml assays.

References:

1. Lowry O. H., et al. (1951) J. Biol. Chem **193**, 265 – 275

DANGER

Hazard

Statements

H314 - H412

Precautions

P260 - P273 - P280

Reaction

P301 +P330 +P331 - P303 +P361 +P353 - P304 +P340 - P305
+P351 +P338 - P312

Storage

P405

HS: 38229000

Storage Temperature: +2 °C to +8 °C