

Product Specification

Catalogue No:	5100-2L
Product Description:	2000ml Conical Flask / Wide Neck / Heavy Duty Rim
Brand:	PYREX
Type of Glass:	Borosilicate glass Type 3.3, which conforms to ISO3585 European and US Pharmacopoeia Type 1glass
Printing:	Product branding, graduation lines and figures and a marking spot (Colour white)
Product Specification:	Overall Height = 282mm Body Diameter = 157mm Neck Internal Diameter to Suit a BS stopper size 58
Pack Quantity:	Quantity 1 pieces per pack

SciLabware Limited confirms that the above mentioned product has been manufactured to meet all relevant Standards as detailed by our Quality Management System.

M Ashley
(Quality Technician)