

## TECHNICAL DATA SHEET

---

<b>Part Nbr</b>	232522
<b>Designation</b>	Pyrex Flask



### Description :

- ✓ Ideal for use in preparation of cultures
- ✓ Compatible with shaker equipment
- ✓ Manufactured from Pyrex borosilicate glass with white enamel markings
- ✓ Compatible with either aluminium or polypropylene closures
- ✓ Capacity : 250 ml
- ✓ Dimensions (body  $\varnothing$  x neck  $\varnothing$  x height) 85 x 38 x 148 mm

### Sales unit :

10 Pyrex Flasks

---