

TECHNICAL DATA SHEET

Part Nbr	232523
Designation	Pyrex Flask



Description :

- ✓ Ideal for use in preparation of cultures
- ✓ Compatible with shaker equipment
- ✓ Manufactured from Pyrex borosilicate glass with white enamel markings
- ✓ Compatible with either aluminium or polypropylene closures
- ✓ Capacity : 500 ml
- ✓ Dimensions (body \varnothing x neck \varnothing x height) 105 x 38 x 183 mm

Sales unit :

5 Pyrex Flasks
