

Product Specification / Certificate of Conformity

Product Code	BNH1000P
Product Description	Flacon col large HDPE, Azlon®
Closure Material	Polypropylene
Material particular	HPDE
Colour	Naturel
Closure Colour	Naturel
Capacity ml	1000
Pack quantity	5

Dimensions and Tolerances

Diameter mm	90
Neck ID mm	22
Height mm	203

Product Compliance

Latex	Yes
Reach	Yes
TSE_BSE Compliant	Oui
Food Compliant	EU 10/2011 / FDA CFR 21 177

This is to certify that the goods detailed above have been inspected and tested in accordance with our BS EN ISO9001 and BS EN ISO13485 quality systems, and conform in all respects to the acceptable quality levels (AQL) detailed therein and any relevant specifications or drawings, unless recorded below as concessionary goods.



Edel Tuite
Head of Quality

Date Issued: 03/07/2019

