

# Product Specification — CO-RE Tip

## P/N and Product description

---

- P/N 235541
- 1000µl wide bore (ø 3.2mm) High CO-RE Filter Tip, non-sterile, black, for Microlab STAR Line



## Material

---

- Polypropylen
- Filter tip insert: Polyethylene
- conductive

## Packaging & processing

---

- 96 tips/tray, 5 trays/blister, 3840 tips/case
- Color-code label on tip tray: white
- Barcode label contains: unique identification number, Lot No., production date, cavity no.
- Barcode type: Code 128 subset B or C
- Tips usable on CO-RE single channel and 96 Probe Head



## Certification / Purity grade

---

- ISO Class 8 Clean Room production, compliant with ISO 9001, GMP and FDA
- CE IVD conform



- DNase-/RNase-free
- Pyrogen-free (endotoxin-free)
- PCR inhibitor-free

## Labeling

---

Each box of disposable tips provides a label which indicates a barcode, part number (REF) and Lot No.

## Measurement

---

- Tip length: 80.1mm
- Inner diameter tip orifice: 3.2mm

**HAMILTON** 

Web: [www.hamiltonrobotics.com](http://www.hamiltonrobotics.com)

Email: [infoservice@hamiltonrobotics.com](mailto:infoservice@hamiltonrobotics.com)

United Kingdom & Ireland  
Tel: +44 (0)121-717-0199

China  
Tel: +86-21-6164-6567

France  
Tel: +33 (0)1 69751616

Italy  
Tel: +39-39-689-33-93

Denmark, Norway,  
Sweden, Finland  
Tel: +46 (0) 8 410 273 73

Germany, Switzerland,  
Austria, Benelux  
Tel: +49 (089) 552649-0