

# Technical Datasheet

---

## Article No. 2356

Dichloromethane for HPLC  
min. 99.9 % CH<sub>2</sub>Cl<sub>2</sub> (stab.)

For laboratory use.

Parameter	Value
Appearance	colourless liquid
Boiling point	39-40 °C
Density	1.3266 g/ml
Mol.-Weight	84.93 g/mol
Stabilized with amylene	30-60 mg/kg
Refractive index (20 °C)	1.422-1.426
Water (KF)	max. 100 mg/kg
Non volatile substances	max. 5 mg/kg
Free acid (as HCl)	max. 5 mg/kg
UV-transmittance at 240 nm	min. 60.0 %
UV-transmittance at 250 nm	min. 92.0 %
UV-transmittance at 255 nm	min. 96.0 %
Colour (Hazen)	max. 10
Filtered through 0.2 µm	