

Product Specification — CO-RE Tip

P/N and Product description

- P/N 235677
- 1000µl wide bore (ø 1.2mm) High CO-RE Filter Tip, sterile, black, for Microlab STAR Line



Material

- Polypropylen
- Filter tip insert: Polyethylene
- conductive

Packaging & processing

- 96 tips/tray, 5 trays/blister, 3840 tips/case
- Color-code label on tip tray: white
- Barcode label contains: unique identification number, Lot No., production date, cavity no.
- Barcode type: Code 128 subset B or C
- Tips usable on CO-RE single channel and 96 Probe Head



Certification / Purity grade

- ISO Class 8 Clean Room production, compliant with ISO 9001, GMP and FDA
- Sterility is validated according to ISO 11137
- CE IVD conform



- DNase-/RNase-free
- Pyrogen-free (endotoxin-free)
- PCR inhibitor-free

Labeling

Each box of disposable tips provides a label which indicates a barcode, part number (REF) and Lot No.

Measurement

- Tip length: 91.8mm
- Inner diameter tip orifice: 1.2mm

HAMILTON 

Web: www.hamiltonrobotics.com

Email: infoservice@hamiltonrobotics.com

United Kingdom & Ireland

Tel: +44 (0)121-717-0199

China

Tel: +86-21-6164-6567

France

Tel: +33 (0)1 69751616

Italy

Tel: +39-39-689-33-93

Denmark, Norway,
Sweden, Finland

Tel: +46 (0) 8 410 273 73

Germany, Switzerland,
Austria, Benelux

Tel: +49 (089) 552649-0