

Product Specification – CO-RE Tip

P/N and Product description

- P/N 235965
- 300µl Standard volume CO-RE Tip, stacked, non-sterile, clear, for Microlab STAR Line




Material

- Polypropylen
- conductive


Packaging & processing

- 96 tips/NTR, 5x4 stacks, 11520 tips/case
- Color-code label on NTR: yellow
- Barcode label contains: unique identification number, Lot No., production date, cavity no.
- Barcode type: Code 128 subset B or C
- Tips usable on CO-RE single channel and 96 Probe Head



Certification / Purity grade

- ISO Class 8 Clean Room production, compliant with ISO 9001, GMP and FDA
- CE IVD conform



- DNase-/RNase-free
- Pyrogen-free (endotoxin-free)
- PCR inhibitor-free

Labeling

Each box of disposable tips provides a label which indicates a barcode, part number (REF) and Lot No.

Measurement

- Tip length: 59.9mm
- Inner diameter tip orifice: 0.4mm

HAMILTON 

Web: www.hamiltonrobotics.com

Email: infoservice@hamiltonrobotics.com

United Kingdom & Ireland

Tel: +44 (0)121-717-0199

China

Tel: +86-21-6164-6567

France

Tel: +33 (01) 69751616

Italy

Tel: +39-39-689-33-93

Denmark, Norway,
Sweden, Finland

Tel: +46 (0) 8 410 273 73

Germany, Switzerland,
Austria, Benelux

Tel: +49 (089) 552649-0