

# Technical Datasheet

---

**Article no.: 2384**

**Dichlormethane anhydrous p. a. (min. 99.95 %, stabilized with amylene)**

CH<sub>2</sub>Cl<sub>2</sub>

For laboratory use.

Parameter	Value
Refractive index (20 °C)	1.422 – 1.426
Water (KF)	max. 50 mg/kg
Colour (APHA)	max. 10
Free acids (as HCl)	max. 5 mg/kg
Tetrachloromethane (CCl <sub>4</sub> )	max. 20 mg/kg
Chloroform (CHCl <sub>3</sub> )	max. 50 mg/kg
Stabilized with amylene over molecular sieve 4A	20 – 60 mg/kg

Version No. 02/2024-09-05